

ULTRA PORTABLE X-RAY

THE FUTURE OF HANDHELD RADIOGRAPHY



HIGHLIGHTED SPECS OF ERAY SMART 5HS



FEATURES THAT EMPOWER CARE



HIGH PERFORMANCE IMAGING

Delivers up to 80kW/SmA output for sharp, high-resolution diagnostic images



ULTRA-LIGHTWEIGHT & PORTABLE

Weighing less than 4 kg designed for effortless mobility and easy operation.



SMART SAFETY SYSTEM

Integrated shielding and intelligent exposure control for safer Imaging practices.



PRECISION DOSE MONITORING

Integrated DAP Meter helps monitor and manage radiation exposure efficiently



WIRELESS OPERATION

Cord-free design for true mobility and flexible clinical workflow.



FAST WORKFLOW

Instant image preview and seamless integration with DR systems.



MULTI-LANGUAGE INTERFACE

Supports English, Hind and Tamil for user convenience



REMOTE CONTROL ENABLED

Supports remote operation for enhanced safety and workflow flexibility

ERAY*Gold* FLAT PANEL DETECTOR & ERAY IMAGESUITE ACQUISITION SOFTWARE

Edusoft Healthcare proudly launches ERAY ImageSuite, our proprietary image acquisition software designed to increase productivity by 1.5x, enhance patient care, and eliminate downtime while securing patient data. With features like Auto stitching, DICOM printing, hardware integration, and advanced chiropractic tools, ERAY ImageSuite redefines efficiency and functionality.

For the first time in India, the software integrates CXR TB Screening AI, Bone Fracture AI, and Bone Suppression AI, (as an optional features) making it the most intelligent and versatile imaging solution with unmatched image quality and advanced AI capabilities.




Paired with the ERAY Gold Premium Detector, which offers lightweight design, extended battery support, and seamless connectivity, the ERAY ImageSuite sets a new benchmark in medical imaging excellence.

FLAT PANEL DETECTOR (ERAY GOLD)

ACTIVE DISPLAY AREA	14"x17" / 17"x17"
DIMENSION	383 x 460 x 15mm / 460 x 460 x 15mm
MATRIX	2.5k x 3k / 3k x 3k
SCINTILLATOR	CsI / GOS
GREY SCALE	16 bit
PIXEL	140 µm
INTERNAL STORAGE	Upto 200 images
PANEL WEIGHT	upto 3kg / upto 3.5kg
FACE LOAD/POINT LOAD	300kg/120 kg
IP RATING	IP56
SPATIAL RESOLUTION	3.5lp/mm
DQE at 0.5lp/mm	70%
MTF at 0.5lp/mm	Upto 86%



X-RAY KIT (Including Accessories)



Technical Highlights

MODEL	ERAY SMART 5HS	
X-RAY GENERATOR	POWER OUTPUT	Upto 450W
	Kv RANGE	50 - 80kV (1 step, Variable)
	mA RANGE	2mA - 5mA (1 step, Variable)
	mAs RANGE	0.01 - 10mAs
	EXPOSURE TIME	0.01 sec to 2 sec
COLLIMATOR	SID	43cm x 43cm @100cm SID
	LIGHT SOURCE	LED Lamp
	BEAM LIMITING & CENTERING ARRANGEMENT	Available
X-RAY TUBE	FOCAL SPOT	0.8 mm
	TARGET ANGLE	16 degrees
	ANODE HEAT STORAGE CAPACITY	8000 J /(10000 HU)
	INHERENT FILTRATION	2.5mm Al Eq
POWER REQUIREMENT	LI-POLYMER BATTERY	22.2VDC, 2,250mAh
	NUMBER OF EXPOSURES	About 150 Times (70 kV, 5mAs)
	AC POWER ADAPTER	AC 100-240V, 50/60Hz (25.2V / 1.5A)
	CHARGING TIME OF BATTERY	Upto 4 hours
WEIGHT	MAIN BODY	Less than 4kg
SIZE	MAIN BODY	249mm x 194mm x 194mm
DISPLAY SCREEN(7Inch)	SCREEN-APR BASED EXPOSURE SETTINGS	Available
	DISPLAY OF kV, mA, mAs, EXPOSURE TIME & BATTERY CHARGER INDICATOR	Available
	ERROR MESSAGES	Available
	LANGUAGES	English, Hindi & Tamil
	DAP	Available
	LOW BATTERY CUTOFF FUNCTION	Available
WARRANTY	1 Year	

Edusoft Healthcare Limited

 B 83, Mangolpuri Industrial Area, Phase 2, Mangolpuri, Delhi, 110083

 info@edusofthealth.com | shammi_gupta@edusofthealth.com

 Toll-Free : 1800 120 280 280

 DOMESTIC: +91 8851248073, +91 9870227422 | INTERNATIONAL : +91 8826697816, +91 9810081174

• Noida • Lucknow • Chandigarh • Pune • Bhopal • Chennai • Bengaluru • Kolkata • Hyderabad
• Jaipur • Guwahati • Ahmedabad • Mumbai • Nepal • USA

Copyright © 2026 Edusoft Healthcare Ltd. All rights reserved. Specifications and designs are subject to change without notice. Non-metric weights and measurements are approximate. All data were deemed correct at time of creation. Edusoft is not liable for errors or omissions. All brand, product, service names and logos are trademarks and/or registered trademarks of their respective owners and are hereby recognized and acknowledged.



Scan Code or visit
www.edusofthealth.com
to learn more